

PulsePalOutput update for Open Ephys GUI

Alessio Paolo Buccino¹: alessiob@ifi.uio.no

December 5, 2017

¹ Center for Integrated Neuroplasticity (CINPLA), Department of Biosciences, University of Oslo, Oslo, Norway

New features

1. PulsePalOutput is updated to the new event channel structure.
2. PulsePalOutput allows the user to directly communicate with the Pulse Pal and update all parameters from the GUI.

1 Event channel information

The *Trigger* and *Gate* combo boxes now are updated with information about incoming TTL channels. Input TTL sources are stored in the PulsePalOutput processor and updated with the `updateSettings()` call. In Figure 1 is an example using multiple *Phase Detector* modules. When only one detector is applied the list of incoming inputs contains 8 TTL channels from *Phase Detector 1*. If a new detector module is added, 8 extra channels are added to the lists (*Phase Detector 2*). This feature eases and widens the choice of incoming TTL events to trigger/gate Pulse Pal channels.

Limitations and possible improvements If one wants to trigger a Pulse Pal channel in response to multiple TTL events, it is not possible with the current implementation (built upon previous implementation). A possible solution could be using list boxes (<https://juce.com/doc/classListBox>) instead of combo boxes.




Figure 1: (left) When one detector module is instantiated 8 TTL channels from *Phase Detector 1* are in the combo box. (right) If we instantiate a second detector module, now 16 TTL channels are available, 8 from *Phase Detector 1* and 8 from *Phase Detector 2*.

2 Pulse Pal Parameters settings

The plugin is now a visualizer plugin and the user can set all Pulse Pal parameters from the visualizer (Figure 2).

For each of the 4 channels, one can set: *phase1 duration*, *phase1 voltage*, *phase2 duration*, *phase2 voltage*, *inter phase duration*, *resting voltage*, *inter pulse duration*, *burst duration*, *inter burst interval*, *train duration*, and *train delay* using editable labels.

burst duration, *inter burst interval* appear only when the BURST button is toggled and *phase2 duration*, *phase2 voltage*, *inter phase duration* only when the BIPHASIC button is toggled (e.g. Channel 2 in figure).

LINK TRIG 1 and LINK TRIG 2 toggle buttons link the corresponding channel to the digital input 1 and 2 of the Pulse Pal. The *trigger mode* combo box allows to select the channel trigger mode (0 = normal, 1 = toggle, 2 = pulse gated).

The TTL button sets the channel to TTL parameters (e.g. Channel 4 in figure).

PulsePalOutput checks the consistency of the parameters. Parameters are consistent if the following conditions hold:

- if burst is OFF: the single pulse duration ($\text{phase1} + \text{phase2} + \text{inter phase interval}$ if biphasic is selected) + inter pulse interval < train duration
- If burst is ON: single pulse duration + inter pulse interval < burst duration AND burst duration + inter burst interval < train duration

If parameters are not consistent, they are adjusted in the following way:

- if burst is OFF: train duration = single pulse duration + inter pulse interval + 1
- if burst is ON: burst duration = single pulse duration + inter pulse interval + 1 AND train duration = burst duration + inter burst interval

Limitations and possible improvements Currently, the user can not set custom trains.



Figure 2: The PulsePalEditor visualizer allows to set all parameters of the Pulse Pal from the GUI.